SURFASOFT®

EN - INSTRUCTIONS FOR USE

DA - BRUGSANVISNING

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EN - INSTRUCTIONS FOR USE

SURFASOFT®

Monofilament non-adherent micro-perforated protective contact layer

Surfasoft® is intended to be used by healthcare professionals only in accordance with these instructions.

PACKAGING

Description	Product reference	Packaging
Surfasoft® 10 cm x 15 cm	1015S	Box, 10 pcs
Surfasoft® 20 cm x 33 cm	2033S	Box, 25 pcs
Surfasoft® 50 cm x 65 cm	5065S	Box, 5 pcs

PRODUCT DESCRIPTION – INTENDED PURPOSE

Surfasoft® is a unique transparent dressing made of monofilament polyamide. The rolled threads create a mesh with a pore size of 90 microns, producing a strong, smooth and transparent material with an open draining surface of 41%.

This unique porous pattern allows exudate to evaporate and drain without drying the wound surface, while topical preparations used to treat burns (post-surgery) can pass through to the skin graft.

Surfasoft® can be used as fixation material for skin grafts and/or skin substitutes, as an interface for negative pressure wound therapy (NPWT) and as a wound contact layer.

INDICATIONS

Transparent fixation dressing for mesh grafts (skin grafts) and/or skin substitutes

Surfasoft® is a strong, flexible, inert, non-adherent dressing that reduces the friction and shear forces when used for fixation of skin grafts.

Surfasoft® ensures a tight skin grafts – wound fit for optimal graft ingrowth. The Surfasoft® dressing fixes mesh grafts and promotes wound drainage, while minimizing bacterial adhesion. When used with Surfasoft®, absorbent wound dressings can be removed painlessly and as often as needed. Surfasoft® allows for easy wound inspection, can remain in situ for several days (or till wound healing) and can be removed without damaging the underlying tissue. With Surfasoft® dressings, skin grafts in difficult areas, e.g. the armpit, neck, head or back, have a take rate as high as 100%.

Interface when applying negative pressure wound therapy

Surfasoft® can be used as a protective contact layer for large soft tissue defects, large abdominal defects, protecting fragile and vital structures such as the gastrointestinal tract, heart, exposed tendons, ligaments, prosthetics and more.

Surfasoft® also limits granulation tissue ingrowth in black polyurethane foam and makes dressing changes easier and less painful for the patient.

Wound contact layer

As a non-adherent interface Surfasoft® can be used to protect vulnerable wound surfaces and promotes undisturbed wound healing by preventing adhesion to the wound bed (and granulation tissue damage) or to the newly formed epithelium. Surfasoft® maintains an optimal moist environment for wound healing. Surfasoft® can be used safely with most absorbent dressings thanks to the inert contact material. Dressing changes are also much easier to perform and are less painful for the patient.

We recommend using Surfasoft® if:

- fixing skin grafts with conventional wound dressings proves to be painful, unsatisfactory and time consuming;
- conventional wound dressings disturb the healing process;
- conventional wound dressings allow for inadequate drainage;
- changing conventional wound dressings is painful and traumatic

PRODUCT APPLICATIONS

- Autografts and allografts
- Donor sites
- NPWT indications
- Partial & full-thickness wounds
- Traumatic and chronic wounds
- Clean closed surgical incisions
- Dermatologic lesions
- Skin substitutes

SIDE EFFECTS

Surfasoft® has been clinically tested for many years. No side effects have been reported.

CONTRAINDICATIONS

No contraindications are known. If there are indications that use of Surfasoft® has nonetheless led to negative effects, either for the patient or the area to be treated, please inform the distributor or contact the manufacturer immediately.

INSTRUCTIONS FOR USE

! Remove the product in its paper wrapper from the sterile pouch.

! Open the paper wrapper and take out the Surfsoft® layer.

1. Saturate Surfsoft® with an isotonic saline solution for improved handling.
2. After the placement and fixation of skin grafts on the wound bed, apply Surfsoft® over the graft in its entirety or cut in strips. Carefully cut large pleats and folds with scissors to ensure good contact. Small folds and creases will not diminish the effectiveness of Surfsoft®.
3. Anchor the Surfsoft® using staples, sutures or acrylic adhesive distributed evenly across the wound surfaces and along the edges.
4. Remove excess blood from under the Surfsoft® without moving the graft on the wound bed by wiping the surface firmly with a damp gauze.
5. Apply a moist absorbent dressing of gauze or a non-woven dressing over the Surfsoft®.
6. Using light pressure, fix this dressing layer using an elastic bandage with open structure.
7. Depending on the type of bacteria in the wound bed, the absorbent dressing layer may be impregnated with various antibacterial agents (e.g. 0.5% silver nitrate, 1% acetic acid, 0.5% chloramide, 1% chlorhexidine, 1% betadine, etc.).
8. Depending on the degree of bacterial contamination, the absorbent dressing should be changed at least once (and if necessary, two or three times) a day. The Surfsoft® layer should be left on the wound.
9. The staples or sutures that anchor the Surfsoft® dressing may be removed after 5 days. Surfsoft® can be removed after 7-14 days, after the tissue has been thoroughly soaked in a saline solution, or by showering the area. Surfsoft® may not remain in situ for longer than 28 days.

PRECAUTIONS, WARNINGS, INSTRUCTIONS AND LIMITATIONS

Precautions:

Ensure Surfsoft® is saturated with isotonic saline solution for better handling.

Carefully cut large pleats and folds with scissors to ensure good contact with the wound bed.

Warnings:

Avoid moving the graft on the wound bed when removing excess blood from under Surfsoft®.

Ensure that the absorbent dressing layer is suitable for the type of bacteria in the wound bed.

Do not exceed 28 days of continuous use.

Surfsoft® is intended as a single-use device.

Do not use Surfsoft® if the packaging is damaged or has already been opened.

Any shipment problems must be reported to the supplier immediately upon delivery.

Instructions:

Anchor Surfасoft® with staples, sutures or acrylic adhesive distributed evenly over the wound surface and edges.

Apply a moist wound dressing of gauze or a non-woven dressing over Surfасoft®.

Use light pressure with an elastic bandage to fix the dressing layer.

Change the dressings down to and including the absorbent dressing at least once a day, or as frequently as needed, depending on the degree of bacterial contamination.

Staples or sutures that anchor the Surfасoft® dressing may be removed after 5 days.

Surfасoft® can be removed after 7-14 days, after the tissue has been thoroughly soaked in a saline solution or by showering the area.

Limitations:

Small folds and creases in Surfасoft® do not affect its effectiveness.

Surfасoft® should remain on the wound even when the absorbent dressing layers are changed.

Surfасoft® may only be applied by an authorized doctor or nurse.

Do not exceed 28 days of continuous use.

Risks related to reusing Surfасoft®:

If a sterile **Surfасoft® wound dressing** is **reused**, the following risks may occur:

1. **Loss of Sterility:** Reusing the product can introduce bacteria, fungi, or other contaminants, increasing the risk of infection.
2. **Reduced Functionality:** Reusing the product can reduce its flexibility and become (too) rigid, hence wound coverage and protection could be compromised.
3. **Compromised Absorption & Moisture Control:** Previously absorbed exudate can alter the dressing's ability to manage moisture effectively, leading to maceration or dryness.
4. **Material Degradation:** Reuse may cause structural breakdown, leading to tearing, fiber shedding, or reduced mechanical integrity.
5. **Inconsistent Performance & Safety Risks:** The dressing may no longer provide the intended protective barrier, potentially delaying wound healing or causing further irritation.

STERILIZATION

Surfасoft® is sterilized by gamma-radiation. Do not use Surfасoft® if the packaging is damaged or has already been opened.

STORAGE, SHELF LIFE AND HANDLING INFORMATION

Surfасoft® should be stored at room temperature in a dust-free area designed for storage of sterile medical devices. Upon delivery of the shipment, ensure that the packaging materials are intact. Any shipment problems must be reported to the supplier immediately

OTHER INFORMATION

Immediately notify the manufacturer and/or local competent authority of any incident following the use of Surfasoft®.


DISPOSAL

Dispose as clinical waste.

DIGITAL INSTRUCTIONS FOR USE

A digital copy of the SURFASOFT® instructions can be found on www.haromed.com.

SYMBOLS AND EXPLANATIONS

Symbol	Explanations
	CE Mark
	Manufacturer
	Use By Date
	Catalogue Number
	Batch/Lot Number
	Sterilized Using Irradiation
	Do Not Reuse
	Do Not Re-Sterilize
	Do Not Use If Packaging Is Damaged

Symbol	Explanations
	Keep Away From Sunlight
	Keep Dry
	Temperature Limit
	Consult Instructions for Use or Consult Electronic Instructions for Use
	Caution
	Product does not contain latex
	Medical Device
	Unique Device Identifier
	Single sterile barrier system

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GEGEVENS FABRIKANT:**HAROMED B.V.**