


Product Specification

Zuigslang 200cm + Y stuk

Description	Pieces	Material
Tube	2	PVC
3-way Connector	1	PVC Rigid granule
Tubing connector	2	Polypropylene
HDPE bag	1	HDPE
Seal Pouch	1	Medical paper/plastic
Box	2	Carton

- Product will be delivered sterile
- Product is sterilized by Ethylene Oxide
- MDD class IIa
- The product is LATEX FREE
- Product contains DEHP
- No animal resources are used
- Materials are tested referring EN ISO 10993

Approved by:



Eelco Buding
QA Manager

