

# Vario 8 and Vario 8 c/i

## MULTI-PURPOSE LOW-VACUUM SUCTION PUMPS



### **Versatile, mobile...**

- One pump for all low-vacuum suction needs
- Suitable for a wide range of applications
- Portable and lightweight
- Battery-powered option

### **...and safe by design**

- Reliable QuatroFlex™ technology
- Easy to use and easy to connect
- Simplicity of design permits efficient routine cleaning
- Safety knob eliminates risk of accidental adjustment

## VERSATILE, MOBILE AND SAFE BY DESIGN

The Vario 8 and Vario 8 c/i low-vacuum suction pumps are great all-rounders for hospitals and outpatient clinics. They offer versatility, reliability and mobility, and are suitable for a wide range of low-vacuum applications.

The suction pumps are portable thanks to their compact design and battery-powered option. Medela's proven Swiss technology offers precise vacuum-level settings and ultra-reliable operation, making Vario 8 and Vario 8 c/i dependable devices for long-term use in hospital.

The simple design allows for efficient routine cleaning. Last but not least, the devices' intuitive handling and safety features are easy to understand and use so that clinical staff can get familiar with it in a short time.



### Key features at a glance



#### The right device for all low-vacuum needs

- With Vario 8 and Vario 8 c/i hospitals have the right device for all low-vacuum applications – wherever they take place.
- Battery-powered option.



#### Ideal for use in hospitals and outpatient clinics

- Thanks to their light weight and their carrying handle, the suction pumps can be transported easily.
- The safety knob for vacuum regulation makes the pumps ideal for bedside use.



#### Reliability thanks to proven technology

- The pumps are powered by QuatroFlex™, a proven and vibration-free piston/cylinder drive unit.
- The membrane vacuum regulator guarantees a rapid vacuum build-up and allows precise regulation.



#### Simplicity of design to enhance safety and hygiene

- The Vario 8 and Vario 8 c/i pumps are easy to use and easy to clean.
- The intuitive handling minimizes the need for instruction and training.
- The safety features include a safety knob to prevent unintended vacuum adjustments.



#### Quiet operation

- Thanks to their low noise emission, Vario 8 and Vario 8 c/i reduce disturbance in the operation room or at the patient's bedside.



Vario 8 and Vario 8 c/i  
low vacuum  
–9 kPa/–68 mmHg



AC 3.5 kg  
AC/DC 4.2 kg  
without jar



ISO 9001  
ISO 13485  
CE (93/42/EEC), IIa  
2 years warranty



8 l/min. (+/- 10%)



h x w x d  
380 x 170 x 285 mm

Vario 8 c/i  
(constant/intermittent low-vacuum pump)  
On-times 8s / 16s / 32s  
Off times 4s / 8s / 16s



IP21

CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.

Medical Vacuum Technology  
for Healthcare Professionals

Please contact us or your local Medela  
representative for details.

 Medela AG  
Lättichstrasse 4b  
6341 Baar, Switzerland  
www.medela.com

Local contact